
10/1987

Gebrauchsanweisung
Instruction Sheet
Mode d'emploi

516 01/04

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This piece of apparatus is suitable for measuring small forces and torques such as occur, for example, when measuring the magnetic field in various types of current-carrying coils, exploring the magnetic fields of permanent magnets, measuring the force between electrically charged spheres (Coulomb's law of electrostatics) or upon an electric charge in the field of a parallel plate capacitor, and when determining the electric power with a wattmeter arrangement.

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1. Description

1. The **torsion balance** designed by Professor Schürholz (516 01) consists of a sturdy metal frame supporting between two torsion wires a holder into which various at-

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Notes

1. The five-digit numbers quoted in parentheses refer to the catalogue numbers of the respective apparatus.
2. The data DK ... refer to the experiment descriptions published as a collection of leaflets in the „LEYBOLD PHYSICS LEAFLETS“. Published to date:

	Cat. No.
1st series (72 sheets)	599 21
2nd series (52 sheets)	599 22
3rd series (49 sheets)	599 23
4th series (50 sheets)	599 24

3. The specifications and illustrations are not binding in every detail for the design of the apparatus. It is our policy always to keep our manufacturing programme right up to date so that it makes full allowance for the latest knowledge acquired in all scientific and technical fields.

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tachments can be inserted. For observing the angular displacements, a concave mirror is provided at the upper end of the holder, which reflects a beam of light onto a distant scale. The twist of the wires can be varied by turning a graduated torsion head at the upper end and by adjusting a simple screw clamp at the lower end.

The torsion balance is so designed that individual parts are easy to replace.

The following parts are supplied under Cat. No. 516 01:

- a) 1 frame of nonferrous material, heavy base plate supporting a rod approx. 3 cm in diameter and carrying the upper plate with the torsion head. Torsion head and base plate are each provided with an electrically insulated terminal clamp for holding the torsion wires. One half of the circumference of the torsion head is graduated. In front of the graduation, a pointer is provided as zero mark. Overall height of frame: approx. 72 cm;
- b) 1 Sensitive system, consisting of a metal bar to be inserted lengthwise between two torsion wires. The bar has two transverse bores at right angles to one another, into which various attachments, to be used for the measurements, may be inserted; moreover, means for mounting the damping vane and counterweight are provided at its lower end, and a concave mirror, $10 \times 20 \text{ mm}^2$ in size, focal length approx. 35 cm, can be fitted on its upper end;
- c) 2 Pairs of torsion wires, approx. 20 cm long;
- d) 1 Damping vessel with adjustable support;
- e) 1 Damping vane with adjustable counterweight;
- f) 1 Bar for dynamic calibration, 24 cm long, dia. 6 mm;

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2. Operation

The torsion balance should be set up in a place free from vibration; a solid laboratory bench will suffice as a rule.

The upper torsion wire should be hooked into the torsion head. From this the bar with the mirror towards the top is suspended, from the bar in turn the 2nd torsion wire, and from its lower end the unscrewed part of the screw clamp. Then this part is inserted again into the socket at the base, using little force, and screwed tight.

The tension of the torsion wires does not affect the restoring torque. It should be so adjusted that they give a high-pitched sound upon being plucked. The bar suspended between the torsion wires will not wobble then.

If damping of the torsion balance is required, the damping vane supplied with it can be inserted into the suspended bar and made to dip into the damping vessel arranged on the adjustable support. Depending upon the desired degree of damping, the vessel may be filled up to various levels with distilled water admixed with means to reduce surface tension or with oil.

In order to compensate for the weight of the damping vane, the appropriate counter-weight must be mounted on the bar. In order to make the centre of gravity of the system lie in the axis of the balance, the vane should be fitted, whilst the stand rod supplied with the balance should be inserted through the upper plate of the frame, and the balance should be laid on the side, so

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that the torsion wires are in a horizontal position. The stand rod fitted on top here serves as a support. The counter-weight should then be mounted on the damping vane in such a manner that its weight is compensated.

For lamp-and-scale reading it is advisable to use lamp house (450 60 and 450 51) and single-lens condenser with sliding aperture (460 17). The lamp should be mounted on a stand (300 02/42) using a Leybold boss (301 01) and set up approx. 30 to 40 cm away from the concave mirror, and the latter should be used to image the slit or a wire fastened across the large aperture on a wall or scale (e. g. 516 04). For convenient evaluation of the results it is advisable to select an integer number of metres (e. g. 3 or 4 m) for the distance between mirror and wall.

3. Calibration

For quantitative evaluation of the measurements the torsion balance must be calibrated. The calibration of the balance by one of the two methods described in the following is easy to perform and is, therefore, recommended.

1. Dynamic calibration

The period of oscillation T of the torsional oscillation of a body of known moment of inertia J is used to determine the restoring torque D , i. e. the torque required to produce an angular displacement of 1 rad ($= 57,30^\circ$) using the equation

$$T^2 = 4 \pi^2 J/D.$$

For this the bar for dynamic calibration supplied with the balance is inserted into the system and the period of oscillation (without damping vane) is measured. The moment of inertia of the bar is

$$J = \frac{ml^2}{12} (= 2,72 \times 10^{-4} \text{ kgm}^2 = 2720 \text{ gcm}^2).$$

In a typical experiment, the period was $T = 6.08$ s; this gave for the restoring torque

$$D = 2.9 \times 10^{-4} \text{ Nm.}$$

Here, the moment of inertia of the system itself was neglected, it being less than 1 % of the moment

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$$D_1 = \frac{1}{2 \cdot 100} D.$$

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$$D_1 = 1,45 \cdot 10^{-6} \text{ Nm.}$$

of inertia of the bar for dynamic calibration inserted into it; therefore, it need not be considered.

2. Static calibration

If, for the calibration of the instrument, one cannot assume the equation of oscillation and the equation for calculating the moment of inertia, then a static calibration can be performed as follows: The balance is laid on one side. The stand rod slid through the upper plate here serves as a support. The small bar for static calibration is inserted into the system and adjusted horizontally by turning the torsion head. In order to fix this position accurately, an adjustable mark (e. g. the pointer of the cursor on vertical scale 311 22) should be brought to exactly level with the one end of the calibrating rod. Now the calibrating weight having a mass of 0.5 g should be inserted into one of the two annular grooves provided at a distance of 50 mm from the axis, and the rod should be returned into the previous, marked position by turning the torsion head. Then one of the torsion wires will be twisted, the other one not. The experiment should then be repeated, inserting the calibrating weight into the other annular groove, to find another value for the twist on the torsion head.

In a typical experiment, the two angular displacements of the torsion head were 101° and 97° , respectively. Thus, in order to twist both wires by a mean angle of 97° , a force of 9.81×10^{-3} N would be required to act at a distance of 50 mm from the axis, which gives for the restoring torque:

$$D = \frac{9.81 \times 10^{-3} \times 0.05 \times 57.3}{97} \\ = 2.90 \times 10^{-4} \text{ Nm.}$$

The values obtained by dynamic and by static calibration agree quite well in the above examples. On considering that the angular displacement is doubled by the mirror, it follows that a 1 cm displacement of the light mark on a scale 1 m away corresponds to a torque of

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For more details concerning the calibration, see LEYBOLD Physics Leaflet DC 531.781; a contained in the 4th Series (50 Sheets), Order No. 599 24.

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4. Experiments

Measurement of the magnetic field of current-carrying coils (516 22) using the torsion balance as magnetometer

Measurement of the magnetic fields of circular conductors (516 23)

Measurement of the magnetic fields of permanent magnets (516 21)

Measurement of Coulomb's force between electrically charged spheres (516 20)

Measurement of the electrical output with the torsion balance, the wattmeter coil (516 30) and the solenoids (516 22)

Measurement of force exerted on an electric charge in the homogeneous electric field of a parallel plate capacitor (544 22) with the probe plate (from 516 20).

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